

Creating Beauty with Confidence.

VeinViewer®

VeinViewer® by Christie, the world's only high definition vein visualization technology, helps improve targeting for cosmetic procedures.

Using near-infrared (NIR) light to detect superficial veins, VeinViewer projects the vein pattern directly on the skin. Easily identify and avoid vessels for cosmetic injections, helping to circumvent potential complications and reduce costs. Improve the patient experience, and set your practice ahead of the rest with the cutting edge technology of VeinViewer.

Aesthetics and Technology

- ✓ Harmless NIR and LED light makes VeinViewer perfect for cosmetic injections.
- ✓ Potentially avoid complications.
- ✓ Potentially reduce down time for patients.
- ✓ Increase patient satisfaction and word of mouth referrals.



CHRISTIE 

What can VeinViewer® do for you?



VeinViewer®

VeinViewer provides:

- ✓ The most accurate image of any vein visualization device on the market.
- ✓ A high definition image that shows veins up to 10 mm deep and blood patterns up to 15 mm.
- ✓ Multiple imaging modes, making it ideal for every patient type.
- ✓ Hospital grade durability.
- ✓ A 5-year standard warranty.
- ✓ Ongoing training available.

Make it the standard in your practice.

1-877-SEE-VEIN

Illuminating Revolutionary Patient Care.™

CHRISTIE+